

TSH

Thyroid Stimulating Hormone (ELISA)

REF: DS177701



Intended use

The DiaSino TSH assay is an enzyme-linked immunosorbent assay (ELISA) for the in vitro quantitative determination of thyroid stimulating hormone (TSH) in human serum. The assay is useful in the diagnosis of thyroid or pituitary disorders. For professional use only. **Not intended for newborn screening.**

Summary

Thyroid-stimulating hormone (TSH, thyrotropin) is a glycoprotein having a molecular weight of approx. 30000 daltons and consisting of two subunits. Measurement of the serum concentration oftentimes thyrotropin (TSH), a glycoprotein with a molecular weight of 28,000 daltons and secreted from the anterior pituitary, is generally regarded as the most sensitive indicator available for the diagnosis of primary and secondary (pituitary) hypothyroidism.^{1,2} TSH measurements are equally useful in differentiating secondary and tertiary (hypothalamic) hypothyroidism from the primary thyroid disease. TSH release from the pituitary is regulated by thyrotropin releasing factor (TRH), which is secreted by the hypothalamus, and by direct action of T4 and triiodothyronine (T3), the thyroid hormones, at the pituitary. Increase levels of T3 and T4 reduces the response of the pituitary to the stimulatory effects of TRH. In secondary and tertiary hypothyroidism, concentrations of T4 are usually low and TSH levels are generally low or normal. Either pituitary TSH deficiency (secondary hypothyroidism) or insufficiency of stimulation of the pituitary by TRH (tertiary hypothyroidism) causes this. The TRH stimulation test differentiates these conditions. In secondary hypothyroidism, TSH response to TRH is blunted while a normal or delayed response is obtained in tertiary hypothyroidism.

Test principle

Sandwich principle. Total duration of assay: **80 minutes**.

- Sample, Anti-TSH coated microwells and enzyme labeled Anti-TSH are combined.
- During the incubation, TSH presents in the sample is allowed to react simultaneously with the two antibodies, resulting in the TSH molecules being sandwiched between the solid phase and enzyme-linked antibodies.
- After washing, a complex is generated between the solid phase, the TSH within the sample and enzyme-linked antibodies by immunological reactions.
- Substrate solution is then added and catalyzed by this complex, resulting in a chromogenic reaction. The resulting chromogenic reaction is measured as absorbance.
- The absorbance is proportional to the amount of TSH in the sample.

Reagents

Materials provided

- **TSH Coated Microplate - symbol** TSH PLATE 8 x 12 strips, 96 wells, pre-coated with mouse monoclonal Anti-TSH.
- **TSH Calibrators - symbols** TSH CAL A-F 6 vials, 1 mL each, ready to use; Concentrations: 0(A), 0.5(B), 2.5(C), 7.5(D), 20(E) and 50(F) µIU/mL.
- **TSH Enzyme Conjugate - symbol** TSH CONJ 1 vial, 11 mL of HRP (horseradish peroxidase) labeled mouse monoclonal Anti-TSH in Tris-NaCl buffer containing BSA (bovine serum albumin). Contains 0.1% ProClin300 preservative.
- **Substrate - symbol** SUBSTRATE 1 vial, 11mL, ready to use, (tetramethylbenzidine) TMB.
- **Stop Solution - symbol** STOP 1 vial, 6.0 mL of 1 mol/L sulfuric acid.
- **Wash Solution Concentrate - symbol** WASH 40X 1 vial, 25 mL (40X concentrated), PBS-Tween wash solution.
- **IFU** - 1 copy.
- **Plate Lid**: 1 piece.

Materials required (but not provided)

- Microplate reader with 450nm and 620nm wavelength absorbent capability.
- Microplate washer.
- Incubator.
- Plate shaker.
- Micropipettes and multichannel micropipettes delivering 50µl with a precision of better than 1.5%.
- Absorbent paper.
- Distilled water

Precautions and warnings

- For in vitro diagnostic use only.
- Exercise the normal precautions required for handling all laboratory reagents.
- Disposal of all waste material should be in accordance with local guidelines.
- Do not use reagents beyond the labeled expiry date.
- Do not mix or use components from kits with different batch codes.
- All the specimen and reaction wastes should be considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.
- The Stop Solution contains sulfuric acid, which can cause severe burns. In the event of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. Exposure to 1.0% sodium hypochlorite for 30 minutes may be necessary to ensure effective decontamination.
- Some reagents contain 0.05% or 0.1% ProClin 300 which may cause sensitization by skin contact, which must therefore be avoided. Reagents and their containers must be

disposed of safely. If swallowed, seek medical advice immediately and show this container or label.

- Substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
- For information on hazardous substances included in the kit please refer to the Materials Safety Data Sheet (MSDS), which is available on request.
- Do not smoke, drink, eat or apply cosmetics in the work area.
- Do not pipette by mouth. Wear protective clothing, disposable gloves and eye/face protection when handling samples and reagents. Wash hands after use.
- If any of the reagents comes into contact with the skin or eyes, wash the area extensively with water.

Incident report

Any suspected serious incidents related to this assay shall be immediately reported to DiaSino, DiaSino's Authorized Representative in the EU, and the national competent authorities of the Member States where the users and/or patients are located.

Storage and stability

- Store at 2-8°C.
- Seal and return unused reagents to 2-8°C, under which conditions the stability will be retained for 2 months, or until the labeled expiry date, whichever is earlier.

Specimen collection and preparation

- Human serum is recommended for this assay.
- Cap and store the samples at 18-25 °C for no more than 8 hours. Stable for 7 days at 2-8°C, and 1 month at -20°C. Freeze only once.
- Do not use heat-inactivated samples.
- Sediments and suspended solids in samples may interfere with the test result which should be removed by centrifugation. Ensure that complete clot formation in serum samples has taken place prior to centrifugation.
- Avoid grossly hemolytic, lipemic or turbid samples.

Calibration

The DiaSino TSH ELISA has been standardized against the 2nd IRP WHO Reference Standard 80/558.

Recalibration is recommended when a new reagent lot is used, or the quality controls are out of specified range.

Quality control

Each laboratory should have assay controls at levels in the low, normal, and elevated range for monitoring assay performance. The controls should be treated as unknowns and values determined in every test procedure performed. The recommended controls requirement for this assay are to purchase trueness control materials separately and test them together with the samples within the same run. The result is valid if the control values fall within the concentration ranges printed on the labels.

Wash solution (40X dilution)

Add deionized water to the 40X concentrated Wash Solution Concentrate.

Dilute 25 mL of Wash Solution Concentrate with 975 mL of deionized water to a final volume of 1000 mL. Stable for 2 months at room temperature.

Test procedure

Ensure the patients' samples, calibrators, and controls are at ambient temperature (18-25 °C) before measurement. Mix all reagents through gently inverting prior to use.

- Use only the number of wells required and format the microplates' wells for each calibrator and sample to be assayed.
- Add **25 µL of calibrators or samples** to each well.
- Add **100 µL of enzyme conjugate** to each well.
- Shake the **microplate** gently for **30 seconds** to mix.
- Cover the plate with a plate lid and incubate at **37°C** for **60 minutes**.
- Discard the contents of the micro plate by decantation or aspiration. If decanting, tap and blot the plate dry with absorbent paper.
- Add **350 µL of wash solution**, decant (tap and blot) or aspirate. Repeat 4 additional times for a total of **5 washes**. An automated microplate strip washer can be used. At the end of washing, invert the plate and tap out any residual wash solution onto absorbent paper.
- Add **100 µL of substrate** to each well.
- Cover and **Incubate at ambient temperature (18-25°C)** in the dark for reaction for **20 minutes**. Do not shake the plate after substrate addition.
- Add **50 µL of stop solution** to each well.
- Shake for **15-20 seconds** to mix the liquid within the wells. It is important to ensure that the blue color changes to yellow completely.
- Read the absorbance of each well at **450 nm** (using 620 to 630 nm as the reference wavelength to minimize well imperfections) in a micro plate reader. The results should be read within **30 minutes** of adding the stop solution.

Calculation

- Record the absorbance obtained from the printout of the microplate reader.
- Calculate the mean absorbance of any duplicate measurements and use the mean for the following calculation.
- Plot the common logarithm of absorbance against concentration in µIU/ml for each calibrator.
- Draw the best-fit curve through the plotted points on linear graph paper. Point-to-Point method is suggested to generate a calibration curve.

TSH

Thyroid Stimulating Hormone (ELISA)

The following data is for demonstration only and cannot be used in place of data generations at the time of assay.

Sample	Value (µIU/mL)	Absorbance
Calibrator A	0	0.015
Calibrator B	0.5	0.108
Calibrator C	2.5	0.433
Calibrator D	7.5	1.106
Calibrator E	20	1.877
Calibrator F	50	2.826
Control 1	1.10	0.206
Control 2	8.76	1.184
Sample	1.31	0.239

Limitations - interference

- The assay is unaffected by icterus (bilirubin < 600 µmol/L or < 35 mg/dL), hemolysis (Hb < 0.559 mmol/L or < 0.9 g/dL), lipemia (Intralipid < 1200 mg/dL), and biotin < 94 nmol/L or < 23 ng/mL.
- Criterion: Recovery within $\pm 10\%$ of initial value.
- For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.
- Serum TSH values may be elevated by pharmacological intervention. Domperidone, amiodazon, iodide, phenobarbital, and phenytoin have been reported to increase TSH levels.
- The presence of autoantibodies may induce high molecular weight complexes (macro-TSH) which may cause unexpected high values of TSH.
- Serum TSH values may be elevated by pharmacological intervention. Domperidone, amiodazon, iodide, phenobarbital, and phenytoin have been reported to increase TSH levels.
- A decrease in thyrotropin values has been reported with the administration of propranolol, methimazol, dopamine and thyroxine.
- Patients who have received mouse monoclonal antibodies for either diagnosis or therapy can develop HAMA (human Anti-mouse antibodies). HAMA can produce either falsely high or falsely low values in immunoassays which use mouse monoclonal antibodies.

Measuring range

0.02-50 µIU/mL (defined by the lower detection limit and the maximum of the master curve). Values below the detection limit are reported as < 0.02 µIU/mL. Values above the measuring range are reported as > 50 µIU/mL.

Lower detection limit

0.02 µIU/mL

The detection limit represents the lowest analyte level that can be distinguished from zero. It is calculated as the value lying two standard deviations above that of the lowest standard (master calibrator, standard 1+2 SD, repeatability study, n=21).

Expected values

0.37-5.10 µIU/mL

These values correspond to the 2.5th and 97.5th percentiles of results obtained from a total of 577 healthy test subjects examined.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using DiaSino reagents, pooled human sera, and controls in a modified protocol (EP5-A) of the CLSI (Clinical and Laboratory Standards Institute): 2 times daily for 20 days (n = 40). The following results were obtained:

Sample	Repeatability*			Intermediate precision	
	Mean µIU/mL	SD µIU/mL	CV %	SD µIU/mL	CV %
Human Serum 1	0.057	0.005	8.77	0.0052	9.12
Human Serum 2	0.44	0.024	5.46	0.03	6.82
Human Serum 3	3.11	0.146	4.69	0.20	6.43
PC Universal 1	1.63	0.088	5.40	0.09	5.52
PC Universal 2	7.86	0.321	4.09	0.42	5.34

*Repeatability = within-run precision

Method comparison

A comparison of the DiaSino TSH assay (y) with the Elecsys TSH assay (x) using 188 clinical samples gave the following correlations:

Linear regression

$$y = 1.003x - 0.629$$

$$r = 0.9780$$

The sample concentrations were between approx. 0 and 38 µIU/mL.

Analytical specificity

For the monoclonal antibodies used, the following cross-reactivities were found: LH 0.041 %, FSH 0.001 %, hGH and hCG no cross-reactivity.

diasino

Functional sensitivity

0.05 µIU/mL

The functional sensitivity is the lowest analyte concentration that can be reproducibly measured with an intermediate precision CV of $\leq 20\%$.













Hook effect

There is no high-dose hook effect at TSH concentrations up to 1000 µIU/mL.

References

- Barker, S.B., "Determination of Protein Bound Iodine."
- Journal Biological Chemistry*, 173, 175, (1984).

Symbols

			
In vitro diagnostic medical device	Temperature limit	Consult instructions for use	Catalog number
			
Batch code	Date of manufacture	Use-by date	Contains sufficient for <n> tests
			
Manufacturer	Do not use if package is damaged and consult instructions for use	European Conformity	Authorized representative in the European Community



diasino

DiaSino Laboratories Co., Ltd

No.68 Jingnansi Road National Eco & Tech Development Area
Zhengzhou 450000, P.R. China. info@diasino.com



CMC Medical Devices & Drugs S.L

C/Horacio Lengo N° 18, CP 29006 Málaga, Spain. Tel: +34 9512 14054