



## **DiaSino<sup>®</sup> IL-6, PCT and CRP ELISA**

*A biomarker panel to support early recognition  
and management of sepsis*



**diasino<sup>®</sup>**

# Sepsis represents a significant medical burden

Sepsis is a potentially life threatening complication of infection, trauma or burn injury. **It remains one of the leading causes of death around the world** and is a primary cause of mortality from infection.<sup>1,2</sup> More than 1 in 1,000 people in developed countries develop sepsis each year, representing a major burden on healthcare.<sup>3</sup>

The first indication of possible sepsis is systemic inflammatory response syndrome (SIRS), a condition defined by the presence of two or more clinical indicators (Figure 1). Sepsis is defined as SIRS with evidence or suspicion of infection and if one or more organs become dysfunctional it is defined as severe sepsis. The patient is said to be in septic shock when sepsis is accompanied by refractory hypotension/hypoperfusion.<sup>4</sup>

**The mortality rate of sepsis is very high** and increases depending on the stage and severity of the sepsis; 24% of patients with SIRS die, increasing to over 50% of patients with septic shock (Figure 2).<sup>2</sup>

**I. Confirmation of infection**  
Diagnosis of an infection on the basis of microbiological evidence or clinical criteria

**II. Systemic inflammatory host response (SIRS) (at least 2 criteria)**

- Fever ( $\geq 38^{\circ}\text{C}$ ) or hypothermia ( $\leq 36^{\circ}\text{C}$ ) confirmed by rectal, intravascular or intravesical measurement
- Tachycardia: heart rate  $\geq 90$  bpm
- Tachypnea (frequency  $\geq 20/\text{min}$ ) or hyperventilation ( $\text{PCO}_2 \leq 4.3 \text{ kPa} / \leq 33 \text{ mmHg}$ )
- Leukocytosis ( $\geq 12,000/\text{mm}^3$ ) or leukopenia ( $\leq 4,000/\text{mm}^3$ ) or  $\geq 10\%$  immature neutrophils in differential blood count

**III. Acute organ dysfunction (at least 1 criterion)**

- Acute encephalopathy: reduced alertness, disorientation, agitation, delirium
- Relative or absolute thrombocytopenia: decrease in platelet counts by more than 30 % within 24 hours or a platelet count of less than  $100,000/\text{mm}^3$ . Thrombocytopenia due to acute hemorrhage or immunological causes must be ruled out
- Arterial hypoxemia:  $\text{PaO}_2 \leq 10 \text{ kPa} (\leq 75 \text{ mmHg})$  while breathing ambient air or a  $\text{PaO}_2/\text{FiO}_2$  ratio of  $\leq 33 \text{ kPa} (\leq 250 \text{ mmHg})$  on oxygen administration. A clinically manifest heart or lung disease must be ruled out as a cause of hypoxemia
- Renal impairment: diuresis of  $\leq 0.5 \text{ mL/kg/h}$  for at least 2 hours despite adequate volume resuscitation and/or an increase in serum creatinine level to  $>$  twice the upper limit of normal (ULN)
- Metabolic acidosis: Base excess of  $\leq -5 \text{ mmol/L}$  or lactate concentration of  $> 1.5 \times \text{ULN}$

**Sepsis: criteria I and II**  
Severe sepsis: criteria I, II and III  
Septic shock: criteria I and II, as well as a systolic arterial blood pressure of  $\leq 90 \text{ mmHg}$  for at least 1 hour, or mean arterial pressure of  $\leq 65 \text{ mmHg}$ , or the necessity of vasopressor administration to maintain a target systolic arterial pressure of  $\geq 90 \text{ mmHg}$  or mean arterial pressure of  $\geq 65 \text{ mmHg}$ . Hypotension persists despite adequate volume resuscitation and cannot be explained by other causes

Figure 1: Clinical indicators of SIRS/sepsis and organ-specific markers (from the Guidelines of the German Sepsis Society).<sup>4</sup>




Figure 2: The four categories of inflammatory response to infection are associated with a significant increase in mortality.<sup>2</sup>

# Timely diagnosis and treatment of sepsis to improve patient survival

## Diagnosis of sepsis - Time matters

Clinical signs suggestive of sepsis are non-specific, making diagnosis of sepsis difficult.<sup>5</sup> Blood cultures remain the gold standard for diagnosis of patients with suspected sepsis, however, they have several limitations:<sup>6, 7</sup>

- Require up to 5 days to obtain the results
- Lack sensitivity
- Easy to contaminate
- Cannot assist in early sepsis management decisions

Timely diagnosis and initiation of effective antibiotic treatment and sepsis management has been shown to improve patient outcomes. Almost 80 % of patients with severe sepsis survive if treatment is initiated within 1 hour of diagnosis. Survival declines rapidly if treatment is delayed further, dropping to only 42 % if initiated 6 hours after diagnosis.<sup>8</sup> In addition, it is desirable to prevent unnecessary antibiotic usage in order to reduce treatment costs, reduce the risk of adverse drug reactions and to slow the spread of antibiotic-resistant bacteria.<sup>9, 10</sup>

## Biomarkers in sepsis

There is considerable unmet medical need in sepsis, and biomarkers may have an important clinical role to play (Figure 3). Biomarkers can indicate the presence of sepsis, differentiate bacterial from viral or fungal infection, differentiate local from systemic infection, stratify severity of sepsis, may help to guide antibiotic therapy, provide prognostic information, evaluate response to therapy, predict septic complications and predict the development of organ dysfunction. However, the exact role of biomarkers is yet to be defined.<sup>5</sup>

## DiaSino clinical diagnostic solution to sepsis

DiaSino IL-6, PCT, in combination with CRP, deliver rapid, reliable information about the patient’s immediate inflammatory status and likelihood of bacterial sepsis, which is important for antimicrobial therapy management (see Figure 4).

**Diagnosis and differential diagnosis of infection**

- Infection vs non-infection<sup>5</sup>
- Infection type (bacterial, viral, fungal)<sup>5</sup> • Systemic or local<sup>5</sup>

**Recognition and stratification of severity of sepsis<sup>5</sup>**

**Early and appropriate decision regarding antibiotic therapy decisions**

- Initiation, and discontinuation of antibiotics<sup>9, 11</sup>

Figure 3: Unmet medical needs in the management of sepsis and suspected sepsis.

| Acute inflammatory episode | Clinical indication of sepsis                                    | Differential diagnosis              | Severe sepsis/shock |
|----------------------------|------------------------------------------------------------------|-------------------------------------|---------------------|
|                            | Suspicion/treatment                                              | Characterization of infection*      | Therapy stewardship |
| IL-6                       | Temperature<br>Heart rate<br>Breathing rate<br>Leukocytes<br>CRP | Blood culture<br>PCT<br>IL-6<br>CRP | PCT                 |

Figure 3: DiaSino clinical diagnostic solution to sepsis.  
\* Rapid identification of sepsis pathogens is possible with LightCycler® SeptiFast Test.

# Interleukin-6 (IL-6)

*A key mediator for inflammation and an early alarm signal of infection*

**IL-6 is an early marker for inflammation in sepsis**

IL-6, a key mediator for inflammation and an early alarm signal of infection that becomes elevated as part of the inflammatory response, has emerged as a valuable biomarker in the management of sepsis.<sup>12</sup>

**IL-6 levels predict development of septic complications**

In a study of 1,032 patients with severe trauma, patients who subsequently developed septic complications had the highest IL-6 levels on day 1 following injury.<sup>13</sup> Similarly, in a study of 50 patients following major surgery, IL-6 levels were correlated with the development of septic complications during the first 5 days following surgery (area under the curve [AUC] 0.82; 95 % CI: 0.66 – 0.98), with a sensitivity of 90 % and selectivity of 58 % . Furthermore, when IL-6 levels and clinical indicators were combined, sensitivity and selectivity increased to 100 % and 79 %, respectively.<sup>14</sup>

**IL-6 levels predict severity of sepsis<sup>15</sup>**

Early peak IL-6 levels correlate significantly with the development of SIRS and sepsis. The degree of elevation in IL-6 levels can be used to differentiate SIRS from severe sepsis and septic shock, with higher IL-6 levels correlating with increased severity (Figure 4).<sup>15</sup>

**IL-6 levels are associated with patient outcome and organ dysfunction**

As a marker for systemic inflammation, high IL-6 levels may be predictive of future organ dysfunction.<sup>12</sup> In addition, continually elevated IL-6 levels have been reported to be predictive of mortality in patients with sepsis.<sup>15</sup>

**DiaSino Interleukin-6 (IL-6) ELISA assay**

| DiaSino® Interleukin-6 (IL-6) test characteristics |                                                                                                                                                                                                                                                                                                                                                     |
|----------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Testing time                                       | 40 minutes                                                                                                                                                                                                                                                                                                                                          |
| Test principle                                     | Sandwich principle                                                                                                                                                                                                                                                                                                                                  |
| Calibrators                                        | 0, 10, 50, 250, 1000, 3000 pg/mL                                                                                                                                                                                                                                                                                                                    |
| Sample material                                    | Serum, Li-heparin, K2-EDTA and K3-EDTA plasma                                                                                                                                                                                                                                                                                                       |
| Sample volume                                      | 50 µL                                                                                                                                                                                                                                                                                                                                               |
| Detection limit                                    | 2 pg/mL                                                                                                                                                                                                                                                                                                                                             |
| Measuring range                                    | 2–3000 pg/mL                                                                                                                                                                                                                                                                                                                                        |
| Traceability                                       | WHO Standard NIBSC 1 <sup>st</sup> IS 89/548                                                                                                                                                                                                                                                                                                        |
| Expected values                                    | ≤ 7 pg/mL<br>A testing study using the DiaSino assay on samples from 177 apparently healthy individuals, a reference range up to 7 pg/mL IL-6 (95th percentile) was determined.<br>Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges. |

# Procalcitonin (PCT)

*A clear signal for antibiotic starts and stops, and when managing critically ill patients with sepsis*

## Physiological role of PCT

PCT is a propeptide of calcitonin, produced by parafollicular cells of the thyroid and neuroendocrine cells of the lung and intestine.<sup>18</sup> PCT levels are raised in response to bacterial endotoxins and proinflammatory stimuli, including IL-6, interleukin 1-b (IL-1b) and tumour necrosis factor- $\alpha$  (TNF- $\alpha$ ).<sup>11</sup> PCT has been reported to be a sensitive marker for bacterial infection, with especially high levels being present following an infection of bacterial origin.<sup>19</sup> Conversely, interferon- $\gamma$  (IFN- $\gamma$ ) blocks the production of PCT, resulting in low PCT levels following viral infection.<sup>11</sup>

## PCT can aid in the diagnosis and stratification of sepsis

Evidence supports a value for PCT in differentiating patients with sepsis caused by bacterial infection from those with SIRS due to a non-infectious cause. By combining PCT with clinical indicators, the accuracy of sepsis diagnosis may be improved further.<sup>16</sup> PCT levels increase as sepsis progresses and severity increases, allowing the differentiation of patients with SIRS, sepsis, severe sepsis or septic shock.<sup>11</sup> Furthermore, PCT levels correlate with increased organ dysfunction in patients with sepsis (as measured by the Sequential Organ Failure Assessment [SOFA] score)<sup>13</sup>, and are increased in patients with poorer outcome.

## PCT levels as a tool for antibiotic therapy guidance

Effective antibiotic treatment is reflected by declining PCT values. Consequently, serial determinations of PCT can be used to monitor the course of systemic bacterial infections and to tailor therapeutic interventions more efficiently.<sup>11</sup>

In multiple studies in various conditions, the use of PCT to guide antibiotic therapy not only reduced the mean number of days of antibiotic treatment, but also reduced the number of days within the intensive care unit (and thereby costs), with no detrimental effect on patient outcome.

## DiaSino Procalcitonin (PCT) ELISA assay

| DiaSino® Procalcitonin (PCT) test characteristics |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|---------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Testing time                                      | 40 minutes                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Test principle                                    | Sandwich principle                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Calibrators                                       | 0, 0.5, 1.0, 2.5, 10, 25 ng/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Sample material                                   | Serum, Li-heparin, K2-EDTA and K3-EDTA plasma.                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Sample volume                                     | 25 $\mu$ L                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Detection limit                                   | 0.02 ng/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Measuring range                                   | 0.02-80 ng/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Traceability                                      | Standardized against Roche Elecsys® PCT                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Expected values                                   | <p>&lt; 0.05 ng/mL</p> <p>A study performed with DiaSino PCT assay using 472 samples from apparently healthy males (237) and females (235) revealed the following normal value: 0.05 ng/mL (95th percentile).</p> <ul style="list-style-type: none"><li>• &lt; 0.5 ng/mL represent a low risk of severe sepsis and/or septic shock</li><li>• 0.5-2.0 ng/mL, observe according to patient's medical history</li><li>• &gt; 2.0 ng/mL represent a high risk of severe sepsis and/or septic shock</li></ul> |

# C-Reactive Protein (CRP)<sup>17-24</sup>

*The most sensitive of the acute phase reactants during inflammatory processes*

**CRP is the most sensitive of the acute phase reactants during inflammatory processes**

C-reactive protein (CRP) is the classic acute phase protein in inflammatory reactions. It is synthesized by the liver and consists of five identical polypeptide chains that form a five-membered ring having a molecular weight of 105000 daltons. CRP is the most sensitive of the acute phase reactants and its concentration increases rapidly during inflammatory processes. Complexed CRP activates the classical complement pathway. The CRP response frequently precedes clinical symptoms, including fever.

**CRP levels as a tool for diagnosis of major trauma and severe infection (sepsis)**

In normal healthy individuals CRP is a trace protein with a range up to 5 mg/L. After onset of an acute phase response the serum CRP concentration rises rapidly and extensively. The increase begins within 6 to 12 hours and the peak value is reached within 24 to 48 hours. Levels above 100 mg/L are associated with severe stimuli such as major trauma and severe infection (sepsis). CRP response may be less pronounced in patients suffering from liver disease.

**Clinical applications of CRP testing**

CRP assays are used to detect systemic inflammatory processes; to assess treatment of bacterial infections with antibiotics; to detect intrauterine infections with concomitant premature amniorrhexis; to differentiate between active and inactive forms of disease with concurrent infection, e.g. in patients suffering from SLE or Colitisulcerosa; to therapeutically monitor rheumatic disease and assess anti-inflammatory therapy; to determine the presence of post-operative complications at an early stage, such as infected wounds, thrombosis and pneumonia, and to distinguish between infection and bone marrow rejection. Postoperative monitoring of CRP levels of patients can aid in the recognition of unexpected complications (persisting high or increasing levels).

Measuring changes in the concentration of CRP provides useful diagnostic information about how acute and how serious a disease is. It also allows judgements about the disease genesis. Persistence of a high serum CRP concentration is usually a grave prognostic sign which generally indicates the presence of an uncontrolled infection.

**DiaSino C-Reactive (CRP) ELISA assay**


| DiaSino® C-Reactive (CRP) test characteristics |                                                                                                                                                                                         |
|------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Testing time                                   | 40 minutes                                                                                                                                                                              |
| Test principle                                 | Sandwich principle                                                                                                                                                                      |
| Calibrators                                    | 0, 5, 20, 50, 150, 300 mg/L                                                                                                                                                             |
| Sample material                                | Serum, Li-heparin, K2-EDTA and K3-EDTA plasma.                                                                                                                                          |
| Sample volume                                  | 10 µL (diluted patient sample)                                                                                                                                                          |
| Detection limit                                | 0.1 mg/L                                                                                                                                                                                |
| Measuring range                                | 0.1-300 mg/L                                                                                                                                                                            |
| Traceability                                   | Standardized against against IFCC CRM 470                                                                                                                                               |
| Expected values                                | < 5.0 mg/L<br>Expected values may vary with age, sex, diet and geographical location. Each laboratory should determine its own expected values as dictated by good laboratory practice. |

# Use of biomarker panel

*The use of a biomarker panel has the potential to provide greater information than measurement of single marker alone*

IL-6, PCT and CRP are complementary biomarkers that provide different information throughout the course of sepsis (Figures 4 and 5). IL-6 levels increase rapidly following insult and peak at around 2–6 hours in surgical trauma patients.<sup>13</sup> Following insult, IL-6 stimulates expression of acute phase proteins such as CRP and induces PCT production through the immune response cascade.<sup>11,25</sup> PCT levels increase primarily in the presence of a bacterial infection, and raised PCT levels are highly suggestive of bacterial sepsis. CRP is a later marker of inflammation, with levels peaking 12–48 hours after infection.<sup>13</sup>

The combination of PCT and IL-6 has been shown to have greater sensitivity and specificity for the detection of early-onset neonatal sepsis than either marker alone.<sup>26,27</sup> The use of a biomarker panel has the potential to provide greater information than measurement of single markers alone and should be investigated further.<sup>5,16,28</sup>



**Figure 4:** IL-6, PCT and CRP have distinct kinetic profiles following bacterial infection in patients who have undergone thoracic surgery.

## IL-6

- Early mediator and indicator of inflammatory response, infection and sepsis
- Indicative of inflammation and infection earlier than blood cultures
  - IL-6 levels elevated within 2 hours of insult
- Increased levels correlate with increased sepsis severity and decreased survival

## PCT

- Acute phase protein primarily raised during bacterial infection
- PCT levels can drive treatment algorithm for initiation, modification and discontinuation of antibiotics
  - Elevated within 6 hours of insult, and highly suggestive of bacterial sepsis
- Indicative of infection earlier than blood cultures
- Can differentiate between sepsis and severe sepsis

## CRP

- Late marker for inflammation, infection and sepsis
- Confirms increases in IL-6 and PCT to aid diagnosis of sepsis

**Figure 5:** Complementary role of IL-6, PCT and CRP in the management of sepsis.

## References

1. Slade, E., Tamber, P., & Vincent, J. (2003). The Surviving Sepsis Campaign: raising awareness to reduce mortality. *Crit Care*, 7(1), 1-2.
2. Silva, E., de Almeida Pedro, M., Beltrami Sogayar, A., Mohovic, T., de Oliveira Silva, C., et al. (2004). Brazilian Sepsis Epidemiological Study (BASES study). *Crit Care*, 8(4), R251-260.
3. Issrah Jawad, et al (2012). "Assessing available information on the burden of sepsis: global estimates of incidence, prevalence and mortality." *J Glob Health*. Jun; 2(1): 010404. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3484761/>
4. ACCP/SCCM. (1992). Definitions for sepsis and organ failure and guidelines for the use of innovative therapies in sepsis. *Crit Care Med*, 20(6), 864-874.
5. Pierrakos, C., & Vincent, J. (2010). Sepsis biomarkers: a review. *Crit Care*, 14(R15).
6. Dark, P., Wilson, C., Blackwood, B., McAuley, D., Perkins, G., et al. (2012). Accuracy of LightCycler SeptiFast for the detection and identification of pathogens in the blood of patients with suspected sepsis: a systematic review protocol. *BMJ Open*, e000392.
7. Schuetz, P., Christ-Crain, M., & Mueller, B. (2008). Procalcitonin and other biomarkers for the assessment of disease severity and guidance of treatment in bacterial infections. *Adv Sepsis*, 6(3), 82-89.
8. Kumar, A., Roberts, D., Wood, K., Light, B., Parrillo, J., et al. (2006). Duration of hypotension before initiation of effective antimicrobial therapy is the critical determinant of survival in human septic shock. *Crit Care Med*, 34(6), 1589-1596.
9. Nobre, V., Harbarth, S., Graf, J., Rohner, P. & Pugin, J. (2008). Use of procalcitonin to shorten antibiotic treatment duration in septic patients. *Am J Respir Crit Care Med*, 177, 498-505.
10. Schuetz, P., Briel, M., Christ-Crain, M., Stolz, D., Bouadma, L., et al. (2012). Procalcitonin to guide initiation and duration of antibiotic treatment in acute respiratory infections: an individual patient data meta-analysis. *Clin Infect Dis*, 55(5), 651-662.
11. Schuetz, P., Albrich, W. & Mueller, B. (2011). Procalcitonin for diagnosis of infection and guide to antibiotic decisions: past, present and future. *BMC Med*, 9(107).
12. Ventetuolo, C., & Levy, M. (2008). Biomarkers: Diagnosis and risk assessment in sepsis. *Clin Chest Med*, 29, 591-603.
13. Billeter, A., Turina, M., Seifert, B., Mica, L., Stocker, R., et al. (2009). Early serum procalcitonin, interleukin-6, and 24-hour lactate clearance: useful indicators of septic infections in severely traumatized patients. *World J Surg*, 33(3), 558-566.
14. Mokart, D., Merlin, M., Sannini, A., Brun, J., Delperio, J., et al. (2005). Procalcitonin, interleukin-6 and systemic inflammatory response syndrome (SIRS): early markers of postoperative sepsis after major surgery. *Br J Anaesth*, 34(6), 767-773.
15. Harbarth, S., Holeckova, K., Froidevaux, C., Pittet, D., Ricou, B., et al. (2001). Diagnostic value of procalcitonin, interleukin-6, and interleukin-8 in critically ill patients admitted with suspected sepsis. *Am J Respir Crit Care Med*, 164, 396-402.
16. Harbarth, S., Holeckova, K., Froidevaux, C., Pittet, D., Ricou, B., et al. (2001). Diagnostic value of procalcitonin, interleukin-6, and interleukin-8 in critically ill patients admitted with suspected sepsis. *Am J Respir Crit Care Med*, 164, 396-402.
17. Greiling, H., Gressner, A.M., eds. (1995). *Lehrbuch der Klinischen Chemie und Pathobiochemie*, 3rd ed. Stuttgart/New York: Schattauer Verlag, 234-236.
18. Thomas, L. (2008). *Labor und Diagnose*, 7. Auflage, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, 1010-1021.
19. Burtis, C.A., Ashwood, E.R., eds. (2001). *Tietz Fundamentals of Clinical Chemistry*, 5th ed. Pa: WB Saunders Co, 332-333.
20. Thomas, L., Messenger, M. (1993). *Pathobiochemie und Labordiagnostik der Entzündung*. Lab med 17, 179-194.
21. Young, B., Gleeson, M., Cripps, A.W. (1991). C-reactive protein: A critical review. *Pathology* 23, 118-124.
22. Wasunna, A., Whitelaw, A., Gallimore, R., et al. (1990). C-reactive protein and bacterial infection in preterm infants. *Eur J Pediatr* Mar 149(6):424-427.
23. Vergis, N. (2007). Should CRP be used as a marker of infection in patients with liver cirrhosis? *Clin Lab Int* 6, 12-13.
24. Mackenzie, I., Woodhouse, J. (2006). C-reactive protein concentrations during bacteraemia: a comparison between patients with and without liver dysfunction. *Intensive Care Medicine* 32, 1344-1351.
25. Maruna, P., Frasko, R. & Gürllich, R. (2008). Plasma procalcitonin in patients with ileus. relations to other inflammatory parameters. *Physiol Res*, 57, 481-486.
26. 32 Bender, L., Thaarup, J., Varming, K., Krarup, H., Ellermann-Eriksen, S., et al. (2008). Early and late markers for the detection of early-onset neonatal sepsis. *Dan Med Bull*, 55(4), 219-223.
27. 33 Abdollahi, A., Shoar, S., Nayyeri, F., & Shariat, M. (2012). Diagnostic value of simultaneous measurement of procalcitonin, interleukin-6 and hs-CRP in prediction of early-onset neonatal sepsis. *Mediterr J Hematol Infect Dis*, 4(1), e2012028.
28. 34 Vincent, J., Donadello, K. & Schmit, X. (2011). Biomarkers in the critically ill patient: C-reactive Protein. *Crit Care Clin*, 27, 241-251.

Diasino and DiaSino are used, pending and/or registered trademarks belonging to DiaSino, or one of its subsidiaries or one of its companies. Any other trademark is the property of its respective owner. / DiaSino Laboratories Co., Ltd./ Printed in China, 07-2020.

# diasino

### DiaSino Headquarters

DiaSino Laboratories Co., Ltd.  
No.58, 2nd Avenue  
National Eco & Tech Development Area  
450000 Zhengzhou, China  
Phone: +86 371 5569 5727  
[www.diasino.com](http://www.diasino.com)

### DiaSino USA Office

DiaSino Laboratories, Inc.  
iPOCT Research Center  
228 Park Avenue S #45956  
NY 10003  
New York  
USA

